Catalog No. 00-8402

## **SPOT-LIGHT™ CELL TREATMENT REAGENT**

### **Enzyme Reagent for Cell Preparations prior to ISH Detection**

Lot No.

#### **INTENDED USE**

For use in pretreating cell preparations prior to In Situ Hybridization (ISH) detection.

**QUANTITY:** 40 mL.

#### **REAGENT PROVIDED**

This enzyme reagent contains pepsin and is supplied in an acidic buffer.

#### **STORAGE**

Store at 2-8°C. Do not freeze.

#### **STABILITY**

There are no obvious signs to indicate instability of this reagent. It has been quality controlled to assure consistent and reliable performance. Do not use after the expiration date stamped on container. With proper storage there is no significant loss of performance. If this reagent is stored under any conditions other than those specified, those conditions must be validated by the user.

#### **REAGENT PREPARATION**

Reagent is ready to use.

#### **USAGE**

##### **Cell Preparation**

Slides should be treated with an adhesive before mounting cell preparation. Immerse slides in 2x SSC @ 37°C for 60 minutes.

##### **Enzyme Digestion (Optional)**

Perform enzyme digestion as recommended in the instructions that come with the ISH probe. The following procedure may be used if kit instructions are unavailable.

Cover cell preparation with at least 100 µL of enzyme reagent, incubate in a humidifying chamber for 2 minutes at 37°C. Two minutes is recommended for most cell preparation, however, longer incubation times may be required for some cell types and preparations. Overdigestion can result in lost of nuclei and chromosome structure. Some cell preparations do not require enzyme digestion.

Authorized Representative for IVDD 98/79/EC:

Invitrogen Ltd.  
Inchinnan Business Park  
3 Fountain Drive  
Paisley  
PA4 9RF  
UK

[www.invitrogen.com](http://www.invitrogen.com)

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288 • E-mail: [techsupport@invitrogen.com](mailto:techsupport@invitrogen.com)

PI008402

(Rev 10/08) DCC-08-1089

© 2007 Invitrogen Corporation. All right reserved. These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or [www.invitrogen.com](http://www.invitrogen.com)). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses.

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